

# Group I, Precipitation of Insoluble Chlorides

Thursday, March 15, 2007  
1:57 PM

Purpose: To find which analytical Group I cations are present in an unknown sample

Equations:

- 1) Analytical Group I ( $\text{Ag}^+$ ,  $\text{Pb}^{2+}$ ,  $\text{Hg}_2^{2+}$ )
- 2)  $\text{Hg}_2\text{Cl}(s) \rightarrow \text{Hg}_2^{2+}(aq) + 2\text{Cl}^-(aq)$  Least soluble Chloride
- 3)  $\text{AgCl}(s) \rightarrow \text{Ag}^+(aq) + \text{Cl}^-(aq)$
- 4)  $\text{PbCl}_2(s) \rightarrow \text{Pb}^{2+}(aq) + 2\text{Cl}^-(aq)$  Most soluble chloride
- 5)  $\text{PbCrO}_4(s) \rightarrow \text{Pb}^{2+} + \text{CrO}_4^{2-}$
- 6)  $\text{Hg}_2\text{Cl}_2 + \text{NH}_3 \rightarrow \text{Hg} + \text{HgNH}_2\text{Cl} + \text{NH}_4^+ + \text{Cl}^-$
- 7)  $[\text{Ag}(\text{NH}_3)_2]^+ + \text{Cl}^- + 2\text{H}^+ \rightarrow \text{AgCl}(s) + \text{white precipitate}$

Procedure: Listed from handout "Group I, Precipitation of Insoluble Chlorides" by unknown author.

Hazards: Lead, mercury, and chromate are toxic

Data/Observations:

Unknown #grp, 1-109

A. Precipitation of Group I Chlorides

When added 12M HCl and 6M nitric acid, it fizzed and there was a white precipitate.

B. Separation and Identification of Lead Ion

When added  $\text{K}_2\text{CrO}_4$  to solution a mustard yellow precipitate formed.

C. Separation of Silver from Mercury, and Identification of Mercury (I) Ion

When added 5 additional drops of ammonium hydroxide, no black or gray residue formed.

D. Identification of Silver Ion

After 2 drops of 6M nitric acid, the blue litmus paper turned pink.

No precipitate or cloudiness observed.

Conclusion:

We found that our unknown # 1-109 only contained lead since we observed a yellow precipitate in step A.